

# Xtend® Connect

# Technical Data Specification

Xtend Connect is a patented quick adapter that allows the user to disassemble his/her prosthesis in parts to make it easier to change feet, change trousers or just save space, for example when traveling. Otherwise, Xtend Connect works just like any adapter and has the common standard interfaces to fit most prosthetic products on the market.



### Quick fit units and Base units

Quick III dilles	Item number	Description
	itelli lidilidei	Description
E.	A100-Q30-150	Xtend® Connect quick fit unit with tube clamp 30 mm, 150 kg/330 lb
HIR	A100-Q4-150	Xtend® Connect quick fit unit with 4-hole adapter, 150 kg/330 lb
	A100-QF-150	Xtend® Connect quick fit unit with female pyramid adapter, 150 kg/330 lb
	A100-QP-150	Xtend® Connect quick fit unit for prong adapter, 150 kg/330 lb
	A100-BF-150	Xtend® Connect base unit with female pyramid adapter, 150 kg/330 lb
11	A100-BM-150	Xtend® Connect base unit with male pyramid adapter, 150 kg/330 lb
<b>3</b>	A100-B30-150	Xtend® Connect base unit with 30 mm tube, 150 kg/330 lb NEW!

### Xtend Connect Configurator




Find out what components you need with our digital Connect Selection tool!

Not sure what components that are needed for a specific configuration? Scan the QR-code and use our online configurator.

# PPR0008\_US-04\_2022-10-21 | We reserve the right to make technical changes.

# Technical facts - combinations



**A** Functional height in mm **B** Effective height in mm **C** Total height in mm **D** Effective weight in g **E** Total weight in g \*B30: Adaptable effective height and weight as tube can be cut off with up to 100 mm.

## How to order

The Xtend Connect is delivered with black locking wheel as a standard. For orange locking wheel add -0 as a suffix to the item number. Maximum patient weight is 150 kg/330 lb. Always check that the patient's weight is correct. The maxi- mum weight must never be exceeded. For intensive sports users, the maximum weight is 100 kg/220 lb.

